

## VACOTRON 460

NEW MODEL



The Vacotron 460 is a module for the 4 series. You can apply electrotherapy by using vacuum electrodes with this suction module. This is accomplished fast and simple. Positioning electrodes could not be any easier.

The Vacotron 460 is a valuable addition to the electrotherapy devices of the 4 series.

- Can be connected to all the electrotherapy and combination units of the 4-series
- Quick application of vacuum electrodes on positions where fixation straps cannot be used
- Continuous and intermittent suction
- Accurate adjustment of the suction strength
- Unique cable connection for fast application of rubber electrodes



The Vacotron 460 can only be used in combination with the Sonopuls 492 or the Endomed 482.



# VACOTRON 460

## NEW MODEL

### 1498950 VACOTRON 460 MODULE

Can only be used in combination with the Sonopuls 492 or the Endomed 482



### TECHNICAL SPECIFICATIONS

Dimensions (WxDxH)	24 x 28.6 x 9.3 cm
Weight	2 kg
Vacuum	continuous and pulsed

### STANDARD ACCESSORIES VACOTRON 460

- 1498756 4-Series Information Booklet
- 2x 3444503 Vacuum electrodes Ø 60 mm, set of 2
- 3444505 Sponges Ø 65 mm, set of 4 (for vacuum electrodes Ø 60 mm)
- 2x 3444507 Vacuum lead hose red
- 2x 3444508 Vacuum lead hose black



### VACUUM ACCESSORIES VACOTRON 460

3444509	Vacuum electrodes Ø 30 mm, set of 2
3444503	Vacuum electrodes Ø 60 mm, set of 2
3444504	Vacuum electrodes Ø 90 mm, set of 2
3444516	Sponges Ø 30 mm, set of 4 (for vacuum electrodes Ø 30 mm)
3444505	Sponges Ø 65 mm, set of 4 (for vacuum electrodes Ø 60 mm)
3444506	Sponges Ø 95 mm, set of 4 (for vacuum electrodes Ø 90 mm)
3444507	Vacuum lead hose, red
3444508	Vacuum lead hose, black



### OTHER PRODUCTS THAT MIGHT INTEREST YOU

#### SONOPULS 492

The top device of the 4 series is the Sonopuls 492. This complete unit ensures that you can apply ultrasound therapy, electrotherapy and combination therapy.

#### ENDOMED 482

The Endomed 482 is the electrotherapy device for the demanding therapist.

#### EN-CAR U3

The EN-Car U3 is your companion for any piece of Enraf-Nonius equipment. For laser, ultrasound, biofeedback, shockwave or electrotherapy, with the EN-Car U3 you have it all close at hand.

#### ORDERING DATA

- 1498911 Sonopuls 492 with multi-frequency treatment head, large (5 cm<sup>2</sup>)
- 1498912 Sonopuls 492 with multi-frequency treatment head, small (0,8 cm<sup>2</sup>)
- 1498913 Sonopuls 492 with 2 multi-frequency treatment heads (5 cm<sup>2</sup> en 0,8 cm<sup>2</sup>)

1498920 Endomed 482

1468960 EN-Car U3

Also see the separate product leaflets at [WWW.ENRAF-NONIUS.COM/ELECTROTHERAPY](http://WWW.ENRAF-NONIUS.COM/ELECTROTHERAPY)



**ENRAF-NONIUS**  
 P.O. Box 12080  
 3004 GB Rotterdam  
 The Netherlands  
 T.: +31-(0)10 - 203 06 00  
 E: [info@enraf-nonius.nl](mailto:info@enraf-nonius.nl)  
[www.enraf-nonius.com](http://www.enraf-nonius.com)